OC SERIES OXYGEN CONCENTRATORS

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EN

CS

User Manual

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SYSMED(CHINA) CO.,LTD
11-2-3, No.17 WenSu Street, HunNan New District,
ShenYang, China

GCE[®]
Healthcare

TABLE OF CONTENTS

Glossary.....	2
Special Notes	3
Before Using This Oxygen Concentrator	3
1. Introduction	3
2. Using Your OC Series Oxygen Concentrator	7
3. Cleaning, Maintenance and Servicing	9
4. Alarms and Troubleshooting Guide	9
5. Specifications	11
Appendix A: EMC Information	11
Appendix B: Circuit Diagram	14
Limited Warranty	15

GLOSSARY

	Consult instructions for use		Follow Instruction for Use
	Type BF Applied Part		Temperature limit
	Power on		Power off
	Caution		Fragile, handle with care
	Stacking Limit By Number		Manufacturer
	No open flames		Keep dry
	No smoking		This Way Up
IP21	Drip Proof Equipment		AC Power
	CE certification mark		Suitable for homecare use

SPECIAL NOTES

- Please read this manual carefully before using this product and save it for future reference.
- For assistance with this manual, contact your home care provider.
- Use the flow and duration settings of this product as prescribed by your physician.
- This equipment is not to be used as a life support device and supplies supplemental oxygen only.
- Children and patients with certain disabilities may be unable to understand or communicate product alarms, and therefore must be supervised when using this device.
- It is advised that patients maintain a backup oxygen supply in case of a machine malfunction or power shortage.
- If you experience an adverse reaction when using this device, contact a physician immediately.

EN

BEFORE USING THIS OXYGEN CONCENTRATOR

-  **WARNING:** In high concentrations, oxygen vigorously accelerates combustion. Patients must make all efforts to reduce the risk of fire when undergoing oxygen therapy.
-  **WARNING:** This oxygen concentrator and all its components - including cannulas, connections and power cords - must be kept away from all sources of heat, open flame, sparks and static electricity.
-  **WARNING:** Do not smoke or allow others to smoke while you are undergoing oxygen therapy. Do not smoke or allow others to smoke in the same room as the oxygen concentrator and accessories.
-  **WARNING:** Oils and greases are prone to strong spontaneous combustion when exposed to oxygen under pressure. To avoid the risk of fire and personal injury, do not use oils and greases on or around the oxygen concentrator.
-  **WARNING:** This device is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

1. INTRODUCTION

Your physician has determined that supplemental oxygen is of benefit to you and has prescribed an oxygen concentrator set at a specific flow and duration setting to meet your needs. Do not change the flow or duration settings unless your physician tells you to do so.

Please read and understand this entire manual before using the device.

INTENDED USE

The OC Series Oxygen Concentrator is intended to provide supplemental oxygen to patients requiring oxygen therapy. The device is not intended to be life supporting or life sustaining.

ABOUT YOUR OC SERIES OXYGEN CONCENTRATOR


The OC Series Oxygen Concentrator delivers highly concentrated oxygen to patients requiring oxygen therapy. The device concentrates oxygen from the air for medical uses through a process called pressure swing absorption. Your home care provider will show you how to use your OC Series device and will be able to answer any questions you may have.

PARTS OF YOUR OC SERIES OXYGEN CONCENTRATOR

EN


FIG 1: LCD PANE

EN

FIG 2: BUTTON BOARD

ACCESSORIES AND REPLACEMENT PARTS

Only use approved replacement parts.

The use of incompatible parts or accessories can result in reduced device performance. Contact your home care provider if you have any questions about the use of accessories or replacement parts.

WARNINGS AND CAUTIONS

A warning represents the possibility of harm to the operator or patient.

-  **WARNING:** For your concentrator to operate properly, air must be able to flow freely around the device. The ventilation ports are located at the rear base of the device and at the side air inlet filter. Keep the device at least 30 cm away from walls, furniture and especially curtains that may reduce airflow to the device. Do not place the concentrator in a small closed space (such as a closet).
-  **WARNING:** Do not remove the covers of this device. Servicing must be referred to an authorized and trained GCE home care provider.
-  **WARNING:** In the event of an equipment alarm or if you are experiencing any signs of discomfort consult your home care provider and/or your health care professional immediately.
-  **WARNING:** Oxygen generated by this concentrator is supplemental and should not be considered life supporting or life sustaining. In certain circumstances oxygen therapy can be hazardous; any user should seek medical advice prior to using this device.
-  **WARNING:** It is essential that any patient who is reliant on the oxygen supplied by this device has a backup source of oxygen in place.
-  **WARNING:** Do not use the oxygen concentrator if either the plug or power cord is damaged. Do not use extension cords or electrical adapters.
-  **WARNING:** Do not attempt to clean the concentrator while it is plugged into an electrical outlet.
-  **WARNING:** Do not connect the concentrator in parallel or in series with other oxygen concentrators or other oxygen therapy devices.
-  **WARNING:** Device operation above or outside of the voltage, flow, temperature, humidity and/or altitude values specified may decrease oxygen concentration levels.
-  **WARNING:** Your home care provider is responsible for performing appropriate maintenance at the intervals recommended by the device manufacturer.

CAUTIONS

A caution represents the possibility of damage to the equipment.

- Do not place liquids on or near the device.
- If liquid is spilled on the device, turn the power off and unplug from the electrical outlet before attempting to clean up the spill. Contact your home care provider if the device does not continue to work properly.

DANGER

To reduce the risk of burns, electrocution, fire or personal injury:

- Do not disassemble. There are no user-serviceable parts. Refer all servicing to an approved service professional.
- Avoid using the device while bathing. If your physician has prescribed continuous usage, the concentrator must be located in another at least 3 meters away from the bath.
- Do not come into contact with the concentrator while wet. Do not place or store the device where liquids can spill into the concentrator.
- Do not reach for the product if it has fallen into water. Unplug immediately and contact your product dealer.

RADIO FREQUENCY INTERFERENCE

This equipment has been tested and found to comply with EMC limited specified by IEC/EN 60601-1-2. These limits are designed to provide reasonable protection against electromagnetic interference in a typical medical installation.

TRANSPORT STORAGE AND OPERATING CONDITIONS

	OPERATING	TRANSPORT AND STORAGE
TEMPERATURE	10°C to 37°C (50°F to 98°F)	-30°C to 70°C (86°F to 158°F)
RELATIVE HUMIDITY	20% to 65%, no condensation	15% to 95%, no condensation
ALTITUDE	0 to 1828 meters	N/A
ATMOSPHERIC PRESSURE	80 kPa ~ 101 kPa	80 kPa ~ 101 kPa
ENVIRONMENT	Dry, well-ventilated, dust-free and pollution-free Away from electromagnetic interference	Placed upright and vertical at all times
ELECTRICAL	North America: 115V, 50Hz Rest of the world: 220V, 60Hz	N/A

EN

2. USING YOUR OXYGEN CONCENTRATOR

⚠ WARNING: Do not use extension cords or electrical adapters.

1. Choose a location for your concentrator that allows it to draw in air freely. Make sure the device is at least 30 cm away from walls, furniture, curtains, or any other item that may restrict airflow to the device.
2. Install the concentrator's two stage filter:
 - a. Remove the filter access panel
 - b. Insert the inlet filter and turn clockwise to connect
 - c. Replace the filter access panel
3. Set up the humidifier bottle:
 - a. Press the humidifier connector button to remove the humidifier bottle
 - b. Remove the cap by turning the humidifier bottle clockwise
 - c. Fill the humidifier bottle with distilled or cold boiled water between the min and max markings
 - d. Replace the humidifier cap and tighten until secure
 - e. Plug the humidifier bottle into the humidifier connector to reconnect it to the oxygen concentrator

NOTE: Replace the water in your humidifier bottle on a daily basis.

4. After reading and understanding the contents of this manual, plug the power cord into a grounded wall outlet and press the power switch to the On [I] position.
5. Make sure that the humidifier bottle is secure:
 - a. Use your finger to gently block the oxygen outlet on the humidifier bottle for 20 seconds
 - b. Once the float in the flow meter drops to the bottom of the gauge, remove your finger.
 - c. If the bottle makes a whistling sound - it means that the humidifier is properly secured to your device.
 - d. If you cannot hear this sound, remove the humidifier bottle, unscrew the cap, replace tightly and repeat the above test.
 - e. If no whistling sound is heard after taking these steps, contact your home care provider.
6. Adjust the oxygen flow so that the float mark is centered on the line marking the flow rate prescribed by your physician:
 - a. Turn the flow control dial clockwise to increase the output flow rate.
 - b. Turn the flow control dial anti-clockwise to decrease the output flow rate.

7. Connect your oxygen tube to the oxygen outlet and put on your cannula as directed by your home care provider. If you are using an extension tube, connect this to the oxygen outlet and oxygen tube.
8. When you are not using the oxygen concentrator, press the power switch to the Off [O] position to turn off the power.

EN

OPERATION MODES AND LCD DISPLAY GUIDE

MODE	GUIDANCE	LCD DISPLAY
Start within 5 minutes	In the first five minutes of operation, the oxygen purity indicator, hour, minute and separator flicker at the same time.	 <p>Model OC-S has no oxygen indicator Model OC-E</p>
Continuous operation	The LCD indicates elapsed time and current run time.	
Timed setting	<p>Set the length of your treatment by pressing the up and down buttons on the button board to increase or decrease the duration in 10-minute increments.</p> <p>20 seconds after setting your treatment duration, the timer will automatically begin to count down to the end of your treatment. The oxygen flow will stop once the timer runs out.</p>	
Oxygen parameters display	<p>Following an initial warm-up period of less than 30 minutes, the device should reach the required oxygen concentration level for your flow rate.</p> <p>If oxygen purity is higher than 85%, the LCD will display the high O2 purity bar and concentration percentage.</p>	

3. CLEANING, MAINTENANCE AND SERVICING

Always press the power switch to the Off [O] position and unplug your OC Series Oxygen Concentrator before cleaning.

The exterior of your oxygen concentrator should be cleaned weekly using a damp cloth with a mild household cleaner. Always wipe the exterior dry after cleaning.

Your cannula and humidifier must be cleaned and replaced as instructed by the manufacturer or your home oxygen provider.

EN

You must clean the air filter in the OC Series Oxygen Concentrator at least once a week.

1. Open the access panel to the air filter to the side of the oxygen concentrator
2. Remove the air filters from the device
3. Rinse and allow to dry naturally
4. Once dry, reinstall the air filters

Maintenance and Servicing

The device has no user-serviceable parts, and must only be serviced by an authorized provider.

4. ALARM AND TROUBLESHOOTING GUIDE

The guide below lists alarms and common problems, as well as why they happened and what you can do to resolve them. If you are unable to resolve a problem or if you experience an issue that is not listed below, please contact your home care provider.

ALARMS

ALARM	WHY IT HAPPENED	WHAT TO DO
Malfunction code H:07 is shown on the LCD display.	Device pressure is higher than specification.	Contact your home care provider.
Malfunction code H:06 is shown on the LCD display.	Device pressure is lower than specification.	Contact your home care provider.
Malfunction code H:05 is shown on the LCD display.	The device has overheated.	Contact your home care provider.
Malfunction code E:01 is shown on the LCD display	Output flow is too low.	Contact your home care provider.
Malfunction code E:03 is shown on the LCD display.	OSD sensor circuit failure.	Contact your home care provider.
Power failure alarm.	The device has detected a power outage or electrical supply failure.	Connect to a backup oxygen source and contact your home care provider.

TROUBLESHOOTING GUIDE

The guide below lists common problems, as well as why they happened and what you can do to resolve them. If you are unable to resolve a problem, or if you experience an issue that is not listed below, please contact your home care provider.

EN

PROBLEM	WHY IT HAPPENED	WHAT TO DO
The device is not working when it is turned on, but the LCD display is showing.	Internal part failure.	Contact your home care provider.
The device is not working when it is turned on and the LCD display is not showing or appears intermittently.	The power cord is damaged.	Check whether the power cord is damaged.
	The power cord plug is not properly inserted into the electrical outlet.	Make sure the device is properly plugged in to the electrical outlet.
	The oxygen concentrator is not receiving power from the electrical outlet.	Check your household outlet fuse or circuit.
	Internal part failure.	Contact your home care provider.
Oxygen is not flowing or the oxygen flow is limited.	The oxygen tubing is kinked or blocked, stopping the delivery of oxygen.	Unkink the oxygen tubing or replace tubing if necessary.
	The humidifier bottle is not properly connected to the device.	Reinstall the humidifier bottle.
The float in the flow meter does not move up or down when adjusting the flow control dial.	The flow control dial is not open.	Turn the flow control dial slowly and carefully.
	The flow control dial is faulty.	Contact your home care provider.
There is water in the cannula.	<p>There has been a change in temperature;</p> <p>OR</p> <p>The device is too close to a wall, curtains or furniture.</p>	Dry the inside of the humidifier cover.
		Do not use hot water in the humidifier bottle.
		Do not overfill the humidifier bottle.
		<p>Keep the oxygen concentrator and cannula in the same room at the same temperature.</p> <p>If using an extension hose with your oxygen concentrator, speak to your home care provider about fitting a water trap to collect excess moisture.</p>

5. SPECIFICATIONS

Technical Specifications

	OC-E80	OC-E100
Power consumption (W)	650 750W 6A (North America)	650 750W 6A (North America)
Rating current (A)	2.25	2.25
Electrical requirements	220V (+/- 10%) 50Hz (+/- 10%)	115V (+/- 10%) 60Hz (+/- 10%) (north America)
Flow rate (L/min)	1 to 8	1 to 10
Purity	93±3%	90±3%
Outlet pressure (mpa)	0.05±0.005	0.05±0.005
Dimensions	381 mm x 347 mm x 689 mm	381 mm x 347 mm x 689 mm
Net weight	27 kg	27 kg
Sound level	50 dBA or less	50 dBA or less
Equipment class and type	Class IIa	Class IIa
EMC compliance	IEC 60601-1-2	IEC 60601-1-2
Features	Standard: LCD display; oxygen concentration indicator (red, green, yellow bar on the LCD display); temperature alarm; high/low pressure alarm; low oxygen concentration alarm; power failure alarm; timed setting; intelligent fault diagnosis; maintenance reminder. Additional: Positive pressure outlet	

EN

Disposal

Dispose of this device in accordance with local regulations.

APPENDIX A: EMC INFORMATION

If you have any questions regarding the guidance and declarations listed below, please contact your home care provider.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure that it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT-GUIDANCE
RF emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	
Emission of Radio Frequency Energy RTCA/DO-160G Section 21	Category M	This device is suitable for use onboard commercial airplanes inside the passenger cabin.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT-GUIDANCE
Electrostatic Discharge (ESD) IEC 61000-4-2	±15kV air ±8kV contact	±15kV air ±8kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input-output lines	±2 kV for supply mains ±1 kV for input/output lines	Mains power quality should be that of a typical home or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical home or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical home or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical hospital or home environment.

Note: UT is the a.c. mains voltage prior to application of the test level.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and This Device

Guidance and Manufacturer's Declaration - Electromagnetic Immunity: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure that it is used in such an environment.

IMMUNITY TEST	IEC60601 LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Conducted RF IEC 61000-4-6	3V 150kHz-80MHz 6 Vrms	3V 150kHz-80MHz 6 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d=1.2 \sqrt{P}$ 150kHz-80MHz $d=1.2 \sqrt{P}$ 80MHz-800MHz $d=2.3 \sqrt{P}$ 800MHz-2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	Amateur Radio & ISM Bands between 150kHz and 80 MHz 10V/m 80MHz-2.7GHz	Amateur Radio & ISM Bands between 150kHz and 80 MHz 10V/m 80MHz-2.7GHz	

EN

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from buildings, objects and people.

A: The field strength intensity of stationary transmitters such as wireless cellular/cordless phones and land mobile radios, amateur radio, AM or FM radio and television broadcasts cannot be predicted theoretically with accuracy. To evaluate the electromagnetic environment, an electromagnetic site survey should be considered. If the measured field strength in the device's location is higher than the RF compliance level as stated in the above table, the device should be observed to verify normal operation. If abnormal device performance is observed, additional measures should be taken to relocate or reorient the device.

B: In the frequency range 150 kHz to 80MHz, the field strength should always be lower than 30 v/m.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THIS DEVICE

This device is intended for use in electromagnetic environments in which radiated RF disturbances are controlled. The user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as outlined below, according to the maximum output power of the communications equipment.


RATED MAXIMUM OUTPUT POWER OF TRANSMITTER (WATTS)	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER (METERS)		
	150KHZ~80MHZ $D=1.2\sqrt{P}$	80MHZ-800MHZ $D=1.2\sqrt{P}$	800MHZ-2.5GHZ $D=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters with maximum output power not stated above, the recommended separation distance (d) in meters can be estimated using the formula applicable to the transmitter frequency, where P is the maximum output power rating of the transmitter in watts (W) provided by the transmitter manufacturer.

NOTE 1: At 80 MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from buildings, objects and people.

APPENDIX B: CIRCUIT DIAGRAM

LIMITED WARRANTY

MODEL	WARRANTY PERIOD
OC-E80	2 years
OC-E100	2 years

GCE Healthcare warrants that the system shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for a period of 2 years from the date of sale by GCE Healthcare to the dealer.

EN

Accessories are warranted to be free of defects in materials and workmanship for a period of 90 days from the time of purchase. If the product fails to perform in accordance with the product specifications, GCE Healthcare will repair or replace – at its discretion – the defective material or part. GCE Healthcare will pay customary freight charges from GCE Healthcare to the dealer location only. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to material or workmanship.

GCE Healthcare disclaims all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty is given in lieu of all other express or implied warranties, including the implied warranties of merchantability and fitness for a particular purpose. In addition, in no event shall GCE Healthcare be liable for lost profits, loss of good will, or incidental or consequential damages, even if GCE Healthcare has been advised of the possibility of the same. Some states or provinces do not allow the exclusion or limitation of implied warranties or the disclaimer of incidental and consequential damages. Accordingly, the laws of your state or province may give you additional protections.

To exercise your rights under this warranty, contact your local authorized GCE Healthcare representative.

SYSMED(CHINA) CO.,LTD
11-2-3, No.17 WenSu Street, HunNan
New District, ShenYang, China
Visit the MQ web site at:
<http://www.sysmed.cn/>

Distributed by: GCE s.r.o.
Zizkova 381,
583 01, Chotebor
Czech Republic
www.gcegroup.com

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